4.9 Key Steps in Laboratory Synthesis Answers


- 1. a) For a distillation set up, the condenser is placed horizontally and for heating under reflux, the condenser is placed vertically.
 - b) Reaction mixtures are heated under reflux because it allows prolonged heating without the escape of volatile compounds.
 - c) A sample of the liquid compound could be placed in a round bottom flask with anti-bumping granules and the distillation apparatus set up. When the compound begins to distil its boiling point will be indicated on the thermometer. This can be compared against a literature value (obtained from a reliable source) and if they match then it indicates that the liquid compound is the desired substance.
 - d) The water inlet should be at the bottom of the condenser and the water outlet should be at the top of the condenser.
 - e) Anti-bumping granules are added to ensure that the reaction mixture does not boil too violently.
- 2. Vacuum filtration is much faster than normal filtration and thus saves time in the laboratory.
- 3. a) The solvent must be immiscible with water.

The product must not react with the solvent.

The product must be more soluble in the solvent than in water.

b) A diagram that shows the dichloromethane as the lower layer (as it is more dense).

Remember that the equilibrium constant is given by (product) / (reactant). This indicates that there should be 5 times more solute in the dichloromethane layer than in the water layer.





- 4. a) The ether layer forms the upper layer as it is less dense.
 - b) A partition coefficient of 12, means that the ratio of solute is distributed as 12 parts in the ether layer and 1 part in the water layer. To find the mass of solute in the ether layer:

- c) Many organic compounds are more soluble in ether as it is a non-polar solvent and these compounds are often non-polar (from Higher Chemistry, like dissolves in like).
- 5. a) n = m / GFM

n = 2.36 / 118

n = 0.02 moles

b) NaOH:

 $n = c \times v$

 $n = 1 \times 0.0348 (34.8 \text{cm}^3)$

n = 0.0348 moles

2NaOH: (CH₂COOH)₂

2 : 1

0.0348 moles: 0.0174 moles IN THE AQUEOUS LAYER

At the start there were 0.02 moles, therefore the number of moles in the diethyl ether layer:

0.02 moles - 0.0174 moles = 0.0026 moles

- c) (0.0026 / 0.0174) = 0.149
- d) The extraction process could be carried out twice using two portions of 50cm³ diethyl ether or indeed carrying out four extractions using 25cm³ portions of diethyl ether.

